



**Advanced Wound-Care
&
Personal-Care Products**



ADVANCED WOUND-CARE PRODUCTS

● Post Operative Dressings

Different types of soft conformable dressings for fast and efficient wound healing process. Maximum breathability to reduce skin maceration, unique barrier systems.

● ProFix

Non woven wound dressings

- Maximum Post operative care.
- Provide both comfort and protection.
- Safe and secured dressing fixation over the wound.

Available Sizes

6 cm x 10 cm
 10 cm x 10 cm
 15 cm x 10 cm
 20 cm x 10 cm
 25 cm x 10 cm
 30 cm x 10 cm
 35 cm x 10 cm
 40 cm x 10 cm


● OpiFix

Polyurethane wound dressings

- Combine protection, conformability, and transparency of polyurethane film with the cushioning and absorbency of a non-adherent pad.
- Provide a unique barrier system for an easy application.
- Is a high waterproof dressing.
- Reduce the risk of maceration.

Available Sizes

6 cm x 10 cm
 10 cm x 10 cm
 15 cm x 10 cm
 20 cm x 10 cm
 25 cm x 10 cm
 30 cm x 10 cm
 35 cm x 10 cm
 40 cm x 10 cm



● Drapex

Polyurethane Incise Drapes

- Offering maximum protection, conformability, and transparency during surgical operations.
- Reduce the risk of infection during operations.
- Transparent, thin Anti-glare film drapes, waterproof and complete bacterial barrier.
- Because of the high vapor permeability of the film Drapex stays closely adhered to the skin throughout the operations.
- Stretches and regains to allow skin retraction without trauma to the skin surface.

Available Sizes	10 cm x 15 cm
	15 cm x 28 cm
	40 cm x 42 cm
	45 cm x 55 cm

● Dressing Retention

- Well proven Hypoallergenic adhesive tapes offering maximum strength and rigidity.

● AlphaPlast

Silk (Acetate) Adhesive Tapes

- Manufactured with new technology (Hot melt) avoiding the organic solvents that are used in the conventional method leading to skin allergic reaction.
- Tearable by hand both lengthways and widthways ensuring easy of application.

Available in 6 sizes with range
 From: 1.25 cm x 5 m To 10 cm x 5 m
 From: 1.25 cm x 2 yd To 10 cm x 5 yd

● AlphaStrips

Non woven Pre-cut Adhesive Strips

Combine the maximum performance of dressing retention maximum protection with easy of application.

Available Sizes	3 cm x 10 cm
	6 cm x 10 cm
	2 cm x 15 cm
	4 cm x 15 cm

● AlphaFix

Non woven Adhesive Rolls

Provide Smoothness, comfort, and size flexibility

Available Sizes	18 cm x 5 m
	10 cm x 10 m
	20 cm x 10 m

● PoroPlast

is an adhesive tape made from a transparent micro perforated film coated on one side with hypoallergenic adhesive, which allows tearing lengthway and widthway and vapour permeability.

From: 1.25 cm x 5 yd To 10 cm x 5 yd

● Hospatek

Urine Bags

- Translucent to facilitate checking and calibrating to 2000ml for easy of measuring.
- Leak proof seal.

- Non Return Valve.
- Drainage port at the bottom.
- Individually packed.
- Can be supplied EtO sterilized.



PERSONAL-CARE PRODUCTS

A number of very efficient products for different personal-care applications including Enema Irrigators, Hot / Cold Packs, Pain Relief plasters and Hair Removing strips for Ladies (New).

● Hospatek

Enema Irrigators

- Ideal for personal Hygiene.
- Used for Anal and Vaginal Irrigation.
- Comfortable & Practical.
- Welded by the latest welding technology to ensure a leak proof seal.
- Has a tap to regulate the flow of fluid.
- Available in 2 Colours, Blue and Transparent Bag.



● Hospatek

Hot / Cold Packs

Re-usable Gel bag fo Hot / Cold therapy.

Indications

- Treatment and rehabilitation of muscular, tendon, or ligament injuries.
- Acute and sub-acute injuries, where swelling is present (using cold treatment to creat vasoconstriction).
- Chronic conditions such as muscle hypertonicity or rheumatoid arthritis (treatment with heat).



● TherMax

Capsicum Plasters

Manufactured with a new stripe coating technology for better skin aeration.

Indications

- Rheumatic Pain
- Lumbago.
- Sciatica.
- Neuralgia.



● VERA®

Hair Removing Strips for ladies

An innovative easy, quick, none sensitizing, non painful leaves no residue on the skin, and suitable for all akin types.

Available in 2 forms:

- VERA® Hair Removing Strips: for All Skin Types.
- VERA® Hair Removing Strips: for Sensitive Skin.

