



Orthopedic Products

Afri Medical



These Top Quality Products are made of different types of cloth and nonwoven fabrics both natural and synthetic. All have a smooth finish surface for comfort and safety, with firm support and compression. Available in different sizes to meet all customer needs.



● **BonFix**

Plaster of Paris Bandages

- Made of specially woven cloth, uniformly impregnated with the finest quality, fast setting plaster of paris.
- The bandage is spooled on a round core for superior bandage stability.
- Moisture proof packaging improves the protection of the bandages.
- Provide fast setting, minimum plaster loss, outstanding mouldability, and improved handling.

Available Sizes

- 7.5 cm x 2.7 m x 24 Bandages
- 10 cm x 2.7 m x 24 Bandages
- 15 cm x 2.7 m x 12 Bandages
- 20 cm x 2.7 m x 12 Bandages



New Technology

● **BonPad**

Under Cast Padding

● **BonPad Natural**

- A 100% Viscose soft and absorbent natural orthopedic under cast padding.
- Conform well around body contours, and gives even coverage of the limb.
- Reduce the risk of skin maceration and irritation.

● **BonPad Synthetic**

- Composed of 70% Polyester and 30% Viscose.
- Very smooth for maximum comfort, light in weight, high vapor permeability to minimize skin maceration.

● **BonPad Synthetic Antibacterial**

- The synthetic padding incorporates an antibacterial agent.
- Offers maximum protection against infection, with maximum comfort and conformability.

Available Sizes	5 cm x 2.7 m
	7.5 cm x 2.7 m
	10 cm x 2.7 m
	12.5 cm x 2.7 m
	15 cm x 2.7 m
	20 cm x 2.7 m

● **BonPlast**

Elastic Adhesive Bandages

- Made of high quality cotton and rayon cloth fully spread with stripe adhesive for better breathability.
- Designed to provide firm support and compression.

Indications:

- Varicose veins, ulcers, phlebitis and thrombosis.
- Sprains and strains.
- Sports injuries.
- Retaining dressings.

Available Sizes	5 cm x 4.5 m (Stretched)
	7.5 cm x 4.5 m (Stretched)
	10 cm x 4.5 m (Stretched)

Afri Medical



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com