


Orthopedic Products

Afri Medical



These Top Quality Products are made of different types of cloth and nonwoven fabrics both natural and synthetic. All have a smooth finish surface for comfort and safety, with firm support and compression. Available in different sizes to meet all customer needs.

BonFix
Plaster of Paris Bandages

- Made of specially woven cloth, uniformly impregnated with the finest quality, fast setting plaster of Paris.
- The bandage is spooled on a round core for superior bandage stability.
- Moisture proof packaging improves the protection of the bandages.
- Provide fast setting, minimum plaster loss, outstanding mouldability, and improved handling.

Available Sizes
 7.5 cm x 2.7 m x 24 Bandages
 10 cm x 2.7 m x 24 Bandages
 15 cm x 2.7 m x 12 Bandages
 20 cm x 2.7 m x 12 Bandages

BonCast
Synthetic Casting Bandages

- Made of fiberglass impregnated with polyurethane resin in a shape of rolled bandage.
- Provide smooth finish surface for comfort and safety, very conformable for easy application.
- Light in weight and strong for maximum mobility and durability.
- High moisture vapor permeability, which reduces the risk of skin maceration.
- Recommended in secondary fractures and orthopedic remedial cast, and in cast bracing.
- Outstanding moisture proof packaging.
- Available in different colors.

Available Sizes
 7.5 cm x 3.6 m
 10 cm x 3.6 m
 12.5 cm x 3.6 m



● **BonPad**
Under Cast Padding

● **BonPad Natural**

- A 100% Viscose soft and absorbent natural orthopedic under cast padding.
- Conform well around body contours, and gives even coverage of the limb.
- Reduce the risk of skin maceration and irritation.

● **BonPad Synthetic**


- Composed of 70% Polyester and 30% Viscose.
- Very smooth for maximum comfort, light in weight, high vapor permeability to minimize skin maceration.

● **BonPad Synthetic Antibacterial**

- The synthetic padding incorporates an antibacterial agent.
- Offers maximum protection against infection, with maximum comfort and conformability.

Available Sizes

5 cm x 2.7 m
7.5 cm x 2.7 m
10 cm x 2.7 m
12.5 cm x 2.7 m
15 cm x 2.7 m
20 cm x 2.7 m



New Technology

● **BonPlast**
Elastic Adhesive Bandages

- Made of high quality cotton and rayon cloth fully spread with stripe adhesive for better breathability.
- Designed to provide firm support and compression.

Indications:

- Varicose veins, ulcers, phlebitis and thrombosis.
- Sprains and strains.
- Sports injuries.
- Retaining dressings.

Available Sizes

5 cm x 4.5 m (Stretched)
7.5 cm x 4.5 m (Stretched)
10 cm x 4.5 m (Stretched)

Afri Medical



Head Office: 13, Obour Buildings, Salah Salem St., Cairo, Egypt
P.O. Box 215 Heliopolis, Cairo

Factory: 10th of Ramadan City - Industrial Area C3, Egypt

Tel: +2015 361 684 **Fax:** +2015 364 104

E-mail: info@afrimedical.com

Web Site: www.afrimedical.com